	MedMira Laboratories Inc. Quality Management System			
	Revision/Issue	1/2	Document ID	TT01MSDS
Reveal MedMira HIV Test Tray MSDS	Released By	PB	Page	1 of 3
			Date	10/19/06

IDENTITY	MedMira HIV Test Tray
----------	-----------------------

Section I – Identification of the Substance and Manufacturer

Manufacturer's Name MedMira Laboratories Inc.	Emergency Telephone Number 902-450-1588		
Address 155 Chain Lake Drive, Suite 1	Telephone Number for Information 902-450-1588		
Halifax, Nova Scotia	Date Prepared September 19, 2000	Last Revision Date October 6, 2006	
	Prepared by MedMira Quality Systems Department		
Canada B3S 1B3			

Section II - Hazard Ingredients/Identity Information


Hazardous Components (Specific Chemical Identity; Common Name(s))	OSHA PEL	ACGIH TLV	Other Limits Recommended	%
None	N/A	N/A	N/A	N/A

Section III – Physical/Chemical Characteristics

Boiling Point	Not applicable	Specific Gravity (H ₂ O=1)	Not applicable
Vapor Pressure (mm Hg)	Not applicable	Melting Point	Not available
Vapor Density (AIR=1)	Not applicable	Evaporation Rate (Butyl Acetate=1)	Not applicable
Solubility in Water Non-Soluble			
Appearance and Odor Odorless, white plastic test tray with enclosed white membrane Optional: Blue plastic filter casing with white filter medium (for whole blood tests)			

Section VI – Fire and Explosion Hazard Data

Flash Point (Method Used) Not applicable	Flammable Limits	LEL Not applicable	UEL Not applicable
Extinguishing Media Use an extinguisher appropriate to the surrounding material that is burning.			
Special Fire Fighting Procedures None			
Unusual Fire and Explosion Hazards None			

	MedMira Laboratories Inc. Quality Management System			
	Revision/Issue	1/2	Document ID	TT01MSDS
Reveal MedMira HIV Test Tray MSDS	Released By	PB	Page	2 of 3
			Date	10/19/06

Section V – Reactivity Data


Stability	Unstable	Conditions to Avoid None		
	Stable X			
Incompatibility (Materials to avoid) No information available				
Hazardous Decomposition or Byproducts No information available				
Hazardous Polymerization	May Occur		Conditions to Avoid None	
	Will Not Occur	X		

Section VI – Health Hazard Data

Route(s) of Entry:	Inhalation	Skin X	Ingestion
Health Hazards (Acute and Chronic) None			
Carcinogenicity:	NTP	IARC Monographs	OSHA Regulated
Not applicable			
Signs and Symptoms of Exposure None			
Medical Conditions Generally Aggravated by Exposure None			
Emergency and First Aid Procedures In case of contact, flush the skin with lukewarm water. Wash clothing before reuse. Seek medical attention if irritation persists.			

Section VII – Precautions for Safe Handling and Use

Steps to Be Taken in Case Material is Released or Spilled None required
Waste Disposal Method Follow all federal, state, and local regulations for disposal. MedMira recommended that the tray be autoclaved and/or incinerated after contact with blood products.
Precautions to Be Taken in Handling and Storage None
Other Precautions None

	MedMira Laboratories Inc. Quality Management System			
	Reveal MedMira HIV Test Tray MSDS	Revision/Issue	1/2	Document ID
	Released By	PB	Page	3 of 3
			Date	10/19/06

Section VIII – Control Measures

Respiratory Protection (Specify Type)		
None required		
Ventilation	Local Exhaust	Special
None required		
	Mechanical (General)	Other
Protective Gloves	Eye Protection	
Rubber or plastic gloves	Chemical safety goggles	
Other Protective Clothing or Equipment		
Lab coat		
Work/Hygienic Practices		
Wash thoroughly after handling. MedMira recommends that the test tray be used in a biological safety cabinet due to the use of blood products.		

The statements contained herein are offered for informational purposes only and are based upon technical data that MedMira Laboratories Inc. believes to be accurate. It is intended for use only by persons having the necessary technical skill and at their own discretion and risk. Since conditions and manner of use are outside our control, we make no WARRANTY, EXPRESS OR IMPLIED, OF MERCHANTABILITY, FITNESS OR OTHERWISE.